NADAL® Rota-Adenovirus Test (test cassette)

REF 481915N-10

de	Gebrauchsanweisung	2	cs	Návod k použití	30
en	Instructions for use	6	nl	Gebruiksaanwijzing	34
fr	Instructions d'utilisation	10	ro	Instrucțiuni de utilizare	38
es	Instrucciones de uso	14	bg	Инструкции за употреба	42
it	Istruzioni per l'uso	18	el	Οδηγίες χρήσης	46
pl	Sposób użycia	22	sk	Návod na použitie	51
pt	Instruções de Utilização	26		Symbols	58
				Our Teams	60

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15. Referenzen

- Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
- Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
- Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva, E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

Rev. 0, 2022-05-23 OM

1. Intended Use

The NADAL® Rota-Adenovirus Test is a lateral flow chromatographic immunoassay for the qualitative detection of adenovirus and group A rotavirus antigens in human faecal specimens. The test is intended for use as an aid in the diagnosis of adenovirus and rotavirus infections (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® Rota-Adenovirus Test is intended for professional use only.

2. Introduction and Clinical Significance

Rotavirus is the most common pathogen responsible for acute gastroenteritis, mainly in young children. Its discovery in 1973 and association with infantile gastroenteritis represented a significant advancement in the study of acute non-bacterial gastroenteritis. Rotavirus is transmitted by the faecal-oral route, and has an incubation period of 1-3 days. Although specimens collected between the third and the fifth day of illness are ideal for antigen detection, rotavirus may still be detectable for as long as diarrhoea continues. Rotaviral gastroenteritis may result in mortality for a population at risk, such as infants, the elderly or immunocompromised patients. In moderate climates, rotaviral infections occur mainly in the winter months. Endemics and epidemics affecting thousands of people have been reported. Up to 50% of the analysed specimens from hospitalised children suffering from acute enteric disease were positive for rotavirus. The viruses replicate in the cell nucleus and tend to produce a characteristic cytopathic effect (CPE). Because rotavirus is extremely difficult to culture, it is unusual to isolate the virus for the diagnosis of infection. Instead, a variety of techniques have been developed to detect rotavirus in faecal specimens.

Acute diarrhoeal disease in young children is a major cause of morbidity worldwide and is a leading cause of mortality in developing countries. Research has shown that enteric adenoviruses, such as Ad40 and Ad41, are also one of the primary causes of diarrhoea in children, second only to rotaviruses. These viral pathogens have been identified throughout the world and can cause diarrhoea in children all year round. Infections are most frequently observed in children under two years of age, but have also been found in patients of all ages. Further studies indicate that adenoviruses are associated with 4-15% of all hospitalised cases of viral gastroenteritis.

A rapid and accurate diagnosis of gastroenteritis caused by adenovirus is helpful in establishing the etiology of gastroenteritis and implementing related patient management. Other diagnostic techniques such as electron microscopy (EM) and nucleic acid hybridisation are expensive and labour-intensive. Due to the self-limiting nature of adenoviral infections, such expensive and labour-intensive tests may not be necessary.

3. Test Principle

The NADAL® Rota-Adenovirus Test enables the detection of adenovirus and group A rotavirus antigens through the visual interpretation of colour development on the internal test strip. Specific anti-rotavirus antibodies are immobilised in the test line region 'R' of the membrane. Specific anti-adenovirus antibodies are immobilised in the test line region 'A' of the

membrane. During the test, the specimen reacts with the specific anti-rotavirus and/or anti-adenovirus antibodies which are conjugated to coloured particles and pre-coated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there is a sufficient amount of rotavirus antigens in the specimen, a coloured line will develop in the test line region 'R' of the membrane. If there is a sufficient amount of adenovirus antigens in the specimen, a coloured line will develop in the test line region 'A' of the membrane. The presence of the coloured line indicates a positive result, while its absence indicates a negative result.

The formation of a coloured line in the control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® Rota-Adenovirus test cassettes
- 10 specimen collection tubes containing buffer*
- 10 disposable pipettes for liquid faecal specimens (25 µL)
- 1 instructions for use

*containing the following preservative: ProClin™ 300: <0.03%

5. Additional Materials Required

- Optional: Specimen collection containers
- Optional: Stool collectors
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

Note: All expiry dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the instructions for use prior to testing.
- Do not use the test beyond the expiry date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid the cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.

- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

Notes on specimen collection and storage:

The NADAL® Rota-Adenovirus Test is intended for use with human faecal specimens only.

Faecal specimens should be collected in clean and dry specimen collection containers containing no preservatives or transport media. Urine and the excessive dilution of specimens with toilet water can lead to inaccurate test results.

Viral detection is facilitated by collecting specimens at the onset of diarrhoeic symptoms. It has been reported that the maximum excretion of rotavirus in the faeces of patients with gastroenteritis occurs 3-5 days after the onset of symptoms. With adenovirus, maximum excretion occurs about 3-13 days after the onset of symptoms. If specimens are collected too long after the onset of diarrhoeic symptoms, the antigen quantity may not be sufficient to obtain a positive result, or the antigens detected may not be linked to the diarrhoeic episode.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Faecal specimens can be stored at 2-8°C for up to 48 hours. Diluted faecal specimens in specimen collection tubes can be stored at -20°C for up to 6 months. Bring specimens to room temperature prior to testing.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Specimen preparation

Hold the specimen collection tube upright and remove the applicator stick by unscrewing the blue cap. Be careful not to spill or splash buffer solution from the collection tube.

For solid specimens:

Collect approx. 50 mg of faeces (equivalent to the size of ¼ of a pea) by inserting the applicator stick into at least 6 different areas of the faecal specimen. Avoid scooping the faecal specimen.

Note: Please be sure to insert the applicator stick 6 times in a row into the faecal specimen. Do NOT put the applicator stick back into the tube in between collecting the 6 specimens. Be careful not to spill fluid from the tube. The quality of the specimen and compliance with these instructions may affect your test result.

For liquid specimens:

Holding a pipette vertically, collect a faecal specimen and add 2 drops (approx. 50 µL) to the specimen collection tube containing buffer.

Place the applicator stick, along with the faecal specimen (or without the faecal specimen if the specimen is liquid), back into the specimen collection tube and tightly screw the blue cap closed.

Shake the specimen collection tube vigorously to thoroughly mix the specimen and buffer.

9. Test Procedure

Bring tests, specimens extracted in buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the NADAL® Rota-Adenovirus test cassette from the sealed foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.


2. Place the test cassette on a clean and level surface.

3. Shake the specimen collection tube to ensure good specimen dispersion.

4. Holding the specimen collection tube vertically, unscrew the transparent protective cap and transfer 3 drops of the specimen solution into the specimen well (S) of the test cassette.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

5. Start the timer.

6. Wait for the coloured line(s) to appear. Read the test result after 10 minutes. Do not interpret the result after more than 20 minutes.

10 min

10. Result Interpretation

Positive for rotavirus:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'R'.

Positive for adenovirus:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'A'.

Positive for rotavirus and adenovirus:

A coloured line develops in the control line region 'C'. A coloured line develops in the test line region 'R' and another one develops in the test line region 'A'.

Note: The colour intensity in the test line region 'R' or 'A' may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region 'R' or 'A' should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region 'C'. No lines develop in the test line regions 'R' and 'A'.

Invalid:

The control line 'C' fails to develop. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test cassette:

A coloured line appearing in the control line region 'C' is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Rota-Adenovirus Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of adenovirus and group A rotavirus antigens in human faecal specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of adenovirus and

rotavirus antigens can be determined using this qualitative test.

- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- The NADAL® Rota-Adenovirus Test only detects the presence of adenovirus and rotavirus antigens in specimens and should not be used as the sole criterion for a diagnosis.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of an adenovirus and rotavirus infection.
- In order to prevent pseudo-outbreaks, positive test results in newborns must be confirmed using an alternative testing method (PCR).
- If visible blood is present in faecal specimens, false-positive results cannot be ruled out and must therefore be verified using an alternative test method.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® Rota-Adenovirus Test was evaluated using clinical faecal specimens collected from children and young adults in comparison to another rapid test.

The results are presented in the following tables:

NADAL® Rota-Adenovirus Test (rotavirus)	Another rapid test		
	Positive	Negative	Total
	121	0	121
	0	209	209
Total	121	209	330

Diagnostic sensitivity: >99.9% (96.9% - 100%)*

Diagnostic specificity: >99.9% (98.2% - 100%)*

Overall agreement: >99.9% (98.9% - 100%)*

*95% confidence interval

NADAL® Rota-Adenovirus Test (adenovirus)	Another rapid test		
	Positive	Negative	Total
	103	0	103
	0	227	227
Total	103	227	330

Diagnostic sensitivity: >99.9% (96.4% - 100%)*

Diagnostic specificity: >99.9% (98.3% - 100%)*

Overall agreement: >99.9% (98.9% - 100%)*

*95% confidence interval

Analytical performance

Analytical specificity

Cross-reactivity study

Negative specimens were spiked with the following, potentially cross-reactive organisms and tested at a

concentration of 1.0×10^9 organisms/mL using the NADAL® Rota-Adenovirus Test:

<i>Acinetobacter calcoaceticus</i>	Rotavirus, group C
<i>Acinetobacter spp.</i>	<i>Haemophilus influenzae</i>
<i>Moraxella catarrhalis</i>	<i>Klebsiella pneumoniae</i>
<i>Candida albicans</i>	<i>Neisseria gonorrhoeae</i>
<i>Chlamydia trachomatis</i>	<i>Neisseria meningitidis</i>
<i>E. coli</i>	<i>Proteus mirabilis</i>
<i>Enterococcus faecalis</i>	<i>Proteus vulgaris</i>
<i>Enterococcus faecium</i>	<i>Pseudomonas aeruginosa</i>
<i>Streptococcus</i> , group B	<i>Salmonella choleraesuis</i>
<i>Streptococcus</i> , group C	<i>Staphylococcus aureus</i>
Rotavirus, group B	

No cross-reactivity with the specimens was observed when tested using the NADAL® Rota-Adenovirus Test.

Interference study

Negative and low positive controls spiked with the following potentially interfering substances at the concentrations listed below showed no interference with the NADAL® Rota-Adenovirus Test.

Substance	Concentration
Ascorbic acid	1 mg/mL
Atropine	0.4 mg/mL
Bilirubin	600 mol/L
Haemoglobin	10 mg/mL
Oxalic acid	10 mg/mL
Triglyceride	5 mg/mL

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative as well as low, medium and high rotavirus positive and low, medium and high adenovirus positive specimens using 3 independent lots of the NADAL® Rota-Adenovirus tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 97.4% - 100%). The NADAL® Rota-Adenovirus Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 4 replicates of negative as well as low, medium and high rotavirus positive and low, medium and high adenovirus positive specimens. Testing was performed by 6 operators using 3 independent NADAL® Rota-Adenovirus test lots at 3 different sites on 5 separate days. >99% of the specimens were correctly identified (60/60 correct tests per concentration, 95% confidence interval: 99.1% - 100%). The NADAL® Rota-Adenovirus Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the event of any serious incidents related to the performance of the NADAL® Rota-Adenovirus Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

- Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
- Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.

- Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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15. Bibliographie

- Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
- Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
- Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

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15. Referencias

- Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
- Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
- Cubitt, WD (1982) Rotavirus infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

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15. Bibliografia

- Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
- Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
- Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva. E. D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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15. Referências

- Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
- Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
- Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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15. Reference

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

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15. Referenties

- Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
- Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
- Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva, E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

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15. Referințe

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.

2. Wilhelm I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva, E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

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5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." *Am. J. Clin. Pathol.* 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. *Journal of Microbiology, Immunology and Infection* 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. *American Journal of Infection Control* 2018; 46(12):1411–1413

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- Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- Thomas, Eva, E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

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