NADAL[®] Adenovirus/RSV Combo Rapid Test (test cassette)

REF 491909N-20

de	Gebrauchsanweisung	2	cs	Návod k použití	37
en	Instructions for use	7	fi	Käyttöohje	42
fr	Instructions d'utilisation	12	ro	Instrucțiuni de utilizare	47
es	Instrucciones de uso	17	bg	Инструкции за употреба	52
it	Istruzioni per l'uso	22	el	Οδηγίες χρήσης	57
pl	Sposób użycia	27		Symbols	65
pt	Instruções de Utilização	32		Our Teams	68

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15. Referenzen

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.

2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.

3. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.

4. Hall, C,B,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.

5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.

6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. J.Clin. Microbiol. 257 : 763-767.

7. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.

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1. Intended Use

The NADAL® Adenovirus/RSV Combo Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection and differentiation of adenovirus and respiratory syncytial virus (RSV) antigens in human nasopharyngeal swab specimens. The test is intended for use as an aid in the diagnosis of adenovirus and RSV infections in patients with a suspected RSV or Adenovirus infection showing typical symptoms of a respiratory tract infection (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® Adenovirus/RSV Combo Rapid Test is designed for professional use only.

2. Introduction and Clinical Significance

RSV is the most common respiratory virus in infants and young children. It infects virtually all children by the age of two. In most children, the virus causes symptoms resembling those of the common cold. In infants born prematurely and/or with chronic lung disease, RSV can cause a severe or even life-threatening disease. Prior to the introduction of Synagis®, RSV disease resulted annually in over 125,000 hospitalisations. There was a high mortality risk amongst approximately 2% of infants. Symptoms of RSV can include fever and a runny nose as well as other, more serious symptoms like coughing, difficult and rapid breathing or wheezing.

Similar symptoms may also appear in the case of adenovirus infections. Adenoviruses are pathogens which can cause a variety of different symptoms depending on which tissue is affected. Adenovirus infections are often asymptomatic. Most symptomatic infections occur in children, who frequently exhibit such symptoms as fever, cough and upper respiratory tract infections (such as pharyngitis, otitis media and tonsillitis). In rare cases, adenovirus infections can also lead to severe pneumonia and bronchiolitis in young children. Adenovirus infections are increasingly recognised as a cause of severe respiratory and other clinical illnesses in immunocompromised adults. In addition to the aforementioned symptoms, an adenovirus infection may cause highly contagious keratoconjunctivitis and gastrointestinal tract diseases.

In order to protect vulnerable groups such as infants, young children and immunocompromised individuals and, if necessary, to initiate targeted treatment, it can be useful to differentiate between adenovirus and RSV infections.

3. Test Principle

The NADAL® Adenovirus/RSV Combo Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection and differentiation of adenovirus and RSV antigens in human nasopharyngeal swab specimens.

The NADAL® Adenovirus/RSV Combo Rapid Test contains two internal test strips (left side of the test cassette: RSV Test; right side of the test cassette: Adenovirus Test).

Anti-adenovirus antibodies are immobilised in the test line region 'ADV' of the membrane. Anti-RSV antibodies are immobilised in the test line region 'RSV' of the membrane. A specimen is added to an extraction tube containing buffer in order to release adenovirus or RSV antigens.

During the test, extracted antigens, if present in the specimen, bind to anti-adenovirus antibodies or anti-RSV antibodies conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane chromatographically by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-adenovirus antibodies or by anti-RSV antibodies in the test line regions 'ADV' or 'RSV', respectively. Excess coloured particles are captured in the control line regions 'C'. The presence of a coloured line in the test line region 'ADV' or 'RSV' indicates a positive result. The absence of a coloured line in the test line region 'ADV' or 'RSV' indicates a negative result.

The formation of a coloured line in each control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® Adenovirus/RSV Combo Rapid test cassettes
- Additional material provided according to 93/42/EEC:
Due to possible supply shortages of accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

a) 20 sterile swabs, CE 0197

CITOTEST LABWARE MANUFACTURING Co. Ltd.,
No. 339 Beihai West Road, Haimen, 226100
Jiangsu, P.R. China (authorised EU representative:
WellKang Ltd, Enterprise Hub, NW Business
Complex, 1 Beraghmore Rd., Derry, BT48 8SE,
Northern Ireland)

b) 20 sterile swabs, CE 0197

Hangzhou Yiguoren Biotechnology Co. Ltd., Room
201, Building 5, No. 360, Tianmushanxi Road,
Yuhang Street, Yuhang District, Hangzhou, 311100
Zhejiang, China (authorised EU representative:
BioDomi GmbH, Schokoladenring 13, 65795
Hattersheim am Main, Germany)

- 20 extraction tubes containing buffer (approximately 500 µl each)*
- 20 dropper caps
- 1 reagent holder
- 1 instructions for use

*Buffer contains the following preservative: ProClin™ 300: <0.001%.

5. Additional Materials Required

- Timer

6. Storage & Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze the test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation.

Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

Note: All expiry dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the instructions for use prior to testing.
- Do not use the test beyond the expiry date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid the cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

Nasopharyngeal specimen:

- It is important to obtain as much secretion as possible.
- Insert the swab into the nostril that presents the most secretion under visual inspection, parallel to the palate (not upwards), until resistance is met or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx.
- Gently rub and roll the swab several times.
- Slowly remove the swab while rotating it. Specimens can be collected from both nostrils using the same swab, but it is not necessary to collect specimens from both sides if the tip is saturated with fluid from the first collection.

Note:




- To ensure optimal performance, use the swabs supplied in the test kit.

- Swab specimens should be tested immediately after collection. For best test performance, use freshly collected specimens.
- Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of specimens and lead to inaccurate test results.

9. Test Procedure

Bring tests, specimens, extraction tubes containing buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Hold an extraction tube containing buffer and labelled with the patient or control identification vertically. Peel off the aluminum foil from the tube.
2. Insert the swab with the collected specimen into the tube. Swirl the swab 5 times and press it firmly against the wall and bottom of the extraction tube.
3. Carefully remove the swab whilst compressing the wall of the extraction tube against the swab to squeeze out as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.
4. Attach a dropper cap to the tube.
5. Place the tube into the designated area of the reagent holder.
6. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
7. Place the test cassette on a clean and level surface.
8. Invert the tube and dispense 3 drops (approximately 75 µL) of the extracted solution into each specimen well (S1) and (S2) of the test cassette.
9. Start the timer.
10. Wait for the coloured line(s) to appear. Read the test result after 15 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation

Positive for adenovirus:

Two coloured lines develop in the result area. One line develops in the control line region 'C' and the other line develops in the test line region 'ADV'.

Positive for RSV:


Two coloured lines develop in the result area. One line develops in the control line region 'C' and the other line develops in the test line region 'RSV'.

Note: The colour intensity in the test line region 'ADV' or 'RSV' may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region 'ADV' or 'RSV' should be considered a positive result. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.


Negative:

A coloured line develops in each control line region 'C'. No coloured lines develop in the test line regions 'ADV' and 'RSV'.

**Invalid**

The control line 'C' fails to develop. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

**11. Quality Control**

An internal procedural control is included in the test cassette:

A coloured line appearing in each control line region 'C' is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Adenovirus/RSV Combo Rapid Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection and differentiation of adenovirus and RSV antigens in human nasopharyngeal specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of adenovirus and RSV antigens can be determined with this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- The NADAL® Adenovirus/RSV Combo Rapid Test only detects the presence of adenovirus and RSV antigens in specimens and should not be used as the sole criterion for a diagnosis of an adenovirus or RSV infection.

- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of an adenovirus or RSV infection.
- Positive and negative predictive values are highly dependent on prevalence. The local prevalence should be taken into consideration when interpreting diagnostic test results.
- Positive results do not preclude co-infections with other pathogens.
- Please be very careful when collecting nasopharyngeal swab specimens from children.

13. Performance Characteristics**Clinical performance****Diagnostic sensitivity and specificity**

The NADAL® Adenovirus/RSV Combo Rapid Test was evaluated with clinical nasopharyngeal swab specimens whose status was confirmed using PCR.

The results are presented in the following tables.

		PCR		
NADAL® Adenovirus/ RSV Combo Rapid Test (for adenovirus)		Positive	Negative	Total
	Positive	67	2	69
	Negative	1	214	215
	Total	68	216	284

Diagnostic sensitivity: 98.53% (92.13% - 99.74%)*

Diagnostic specificity: 99.07% (96.69% - 99.75%)*

Overall agreement: 98.94% (96.94% - 99.64%)*

*95% confidence interval

		PCR		
NADAL® Adenovirus/ RSV Combo Rapid Test (for RSV)		Positive	Negative	Total
	Positive	111	1	112
	Negative	2	170	172
	Total	113	171	284

Diagnostic sensitivity: 98.23% (93.78% - 99.51%)*

Diagnostic specificity: 99.42% (96.76% - 99.90%)*

Overall agreement: 98.94% (96.94% - 99.64%)*

*95% confidence interval

Analytical performance**Detection limit**

The 95% detection limit of the NADAL® Adenovirus/RSV Combo Rapid Test for recombinant adenovirus antigen (type 3) is 0.05 ng/mL. The 95% detection limit of the NADAL® Adenovirus/RSV Combo Rapid Test for RSV subtype A strains is 93.6 TCID₅₀/mL and RSV subtype B strains is 11.7 TCID₅₀/mL, and was determined via the dilution of inactivated RSV.

Prozone effect

No prozone effect was observed when testing specimens containing an adenovirus antigen concentration as high as 4.17×10⁵ TCID₅₀/mL or specimens containing an RSV antigen concentration as high as 4.68 × 10⁴ TCID₅₀/mL for RSV subtype A strains and as high as 1.17 × 10⁵ TCID₅₀/mL for RSV subtype B strains.

Analytical specificity

Cross-reactivity study

Adenovirus and RSV negative nasopharyngeal swab specimens spiked with the following potentially cross-reacting substances were tested using the NADAL® Adenovirus/RSV Combo Rapid Test.

Pathogen	Concentration
SARS-CoV-1 Urbani	1×10 ⁶ PFU/mL
SARS-CoV-2 USA-WA1/2020	1×10 ⁶ TCID ₅₀ /mL
Influenza A H1N1 virus	6×10 ⁵ TCID ₅₀ /mL
Influenza A H3N2 virus	9×10 ⁵ TCID ₅₀ /mL
Influenza A H5N1 virus	8×10 ⁴ TCID ₅₀ /mL
Influenza A H7N9 virus	9×10 ⁴ TCID ₅₀ /mL
Influenza B N/A	3×10 ⁵ TCID ₅₀ /mL
Adenovirus type 1	3×10 ⁵ TCID ₅₀ /mL
Adenovirus type 2	2×10 ⁵ TCID ₅₀ /mL
Adenovirus type 3	1×10 ⁵ TCID ₅₀ /mL
Adenovirus type 5	4×10 ⁵ TCID ₅₀ /mL
Adenovirus type 7	1×10 ⁶ TCID ₅₀ /mL
Adenovirus type 55	4×10 ⁵ TCID ₅₀ /mL
Respiratory syncytial virus type A/B	3×10 ⁵ TCID ₅₀ /mL
Human coronavirus 229E	1×10 ⁶ PFU/mL
Human coronavirus OC43	1×10 ⁶ PFU/mL
Human coronavirus NL63	1×10 ⁶ PFU/mL
Human coronavirus HKU1	1×10 ⁶ PFU/mL
MERS-Coronavirus Florida/USA-2/Saudi Arabia.2014	1×10 ⁶ PFU/mL
Parainfluenza virus type 1/2/3/4	1×10 ⁵ TCID ₅₀ /mL
Rhinovirus type A16 N/A	1×10 ⁵ TCID ₅₀ /mL
<i>Legionella pneumophila</i> Bloomington-2	1×10 ⁵ cells/mL
<i>Legionella pneumophila</i> 82A3105	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> K	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> Erdman	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> HN878	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> CDC1551	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> H37Rv	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> 475298 [Maryland (D1) 6B-17]	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> 178 [Poland 23F-16]	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> 262 [CIP 104340]	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> Slovakia 14-10 [29055]	1×10 ⁵ cells/mL
<i>Streptococcus pyogenes</i> typing strain T1	1×10 ⁵ cells/mL
<i>Mycoplasma pneumoniae</i> mutant 22	1×10 ⁵ cells/mL
<i>Mycoplasma pneumoniae</i> FH strain of Eaton Agent	1×10 ⁵ cells/mL
<i>Mycoplasma pneumoniae</i> M129-B7	1×10 ⁵ cells/mL
<i>Chlamydia pneumoniae</i> , AR39	1×10 ⁶ IFU/mL
Human Metapneumovirus 16 type A1, strain IA10-2003	1×10 ⁵ TCID ₅₀ /mL
Enterovirus type 68	4×10 ⁵ TCID ₅₀ /mL
<i>Haemophilus influenzae</i> type b, strain Eagan	1×10 ⁶ CFU/mL
<i>Candida albicans</i> CMCC (F) 98001	1×10 ⁶ CFU/mL
<i>Bordetella pertussis</i> A639	1.4×10 ⁶ CFU/mL
<i>Staphylococcus aureus</i> , strain NCTC 8325	1×10 ⁶ CFU/mL
<i>Staphylococcus epidermidis</i> RP62A (MRSE)	1×10 ⁶ CFU/mL
<i>Pneumocystis jirovecii</i> W303-Pji	1×10 ⁶ CFU/mL
Pooled human nasal wash N/A	14% v/v

No cross-reactivity with the specimens was observed when tested using the NADAL® Adenovirus/RSV Combo Rapid Test.

Interference study

Adenovirus and RSV negative nasopharyngeal swab specimens, nasopharyngeal swab specimens spiked with 10 µL of adenovirus positive control as well as nasopharyngeal swab specimens spiked with 10 µL RSV positive control were spiked with 10 µL of the following potentially interfering substances at the concentrations listed below and showed no interference with the NADAL® Adenovirus/RSV Combo Rapid Test.

	Substance	Concentration
Anti-viral drugs	Zanamivir (influenza)	5 mg/mL
	Oseltamivir (influenza)	10 mg/mL
	Artemether/Lumefantrine (malaria)	50 µM
	Doxycycline hyclate (malaria)	70 µM
	Quinine (malaria)	150 µM
	Lamivudine (antiretroviral medication)	1 mg/mL
	Ribavirin (HCV)	1 mg/mL
	Daclatasvir (HCV)	1 mg/mL
Respiratory specimens	Mucin from bovine submaxillary glands, type I-S	100 µg/mL
	EDTA-anticoagulated human blood	5% (v/v)
	Biotin	100 µg/mL
Nasal sprays or drops	Neosynephrin® (phenylephrine)	10% (v/v)
	Afrin® Nasal Spray (oxymetazoline)	10% (v/v)
	Saline nasal spray	10% (v/v)
Homeopathic allergy relief medicine	homeopathic Zicam®	5% (v/v)
	Allergy Relief Nasal Gel	
	Sodium cromoglycate	20 mg/mL
Anti-inflammatory medication	Olopatadine hydrochloride	10 mg/mL
	Paracetamol (acetaminophen)	199 µM
	Acetylsalicylic acid	3.62 mM
Antibiotics	Ibuprofen	2.425 mM
	Mupirocin	10 mg/mL
	Tobramycin	5 µg/mL
	Erythromycin	81.6 µM
	Ciprofloxacin	30.2 µM

Precision

Repeatability

Repeatability was established by testing 20 replicates of adenovirus negative and positive control specimens as well as 20 replicates of RSV negative and positive control specimens using 3 lots of the NADAL® Adenovirus/RSV Combo Rapid tests at 1 site. >99% of the specimens were correctly identified (20/20 correct tests per control, 95% confidence interval: 83.2% - 100%). The NADAL® Adenovirus/RSV Combo Rapid Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of adenovirus negative and positive control specimens as well as 5 replicates of RSV negative and positive control specimens. Testing was performed by 5 operators using 3 NADAL® Adenovirus/RSV Combo Rapid test lots at 3 different sites on 5 separate days. >99% of the specimens were correctly identified (375/375 correct tests per control, 95% confidence

interval: 99% - 100%). The NADAL® Adenovirus/RSV Combo Rapid Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® Adenovirus/RSV Combo Rapid Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Delewa, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
3. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
4. Hall, C.B.,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.
7. Wadell, G. *Laboratory Diagnosis of Infectious Diseases: Principles and Practices.* New York: Springer-Verlag, Volume II, 1988: 284-300.
8. Infektionen durch Adenovirus - Infektionskrankheiten - MSD Manual Profi-Ausgabe (msdmanuals.com), as of August 8th, 2024
9. Adenoviren: www.infektionsschutz.de, as of August 8th, 2024

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treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.

6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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15. Bibliographie

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.

2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.

3. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.

4. Hall, C,B,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and WW.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.

5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol

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15. Bibliografia

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
3. Hall, C.B, R.G. Douglas, and J.M. Gelman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
4. Hall, C.B.,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and WW.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. J.Clin. Microbiol. 257 : 763-767.
7. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
8. Infektionen durch Adenovirus - Infektionskrankheiten - MSD Manual Profi-Ausgabe (msdmanuals.com), as of August 8th, 2024
9. Adenoviren: www.infektionsschutz.de, as of August 8th, 2024

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15. Referências

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C. Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
3. Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
4. Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S. Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N. Engl. J. Med. 308 : 1443-1447.
5. Taber, L.H.V., Knight, B.E., Gilbert, H.W., McClung, S.Z., Wilson, H.J., Norton, J.M., Thurson, W.H., Gordon, R.L., Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
6. Ahluwalia, G.J., Embree, P., McNicol, B., Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. J. Clin. Microbiol. 257 : 763-767.
7. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
8. Infektionen durch Adenovirus - Infektionskrankheiten - MSD Manual Profi-Ausgabe (msdmanuals.com), as of August 8th, 2024
9. Adenoviren: www.infektionsschutz.de, as of August 8th, 2024

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15. Reference

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
3. Hall, C.B. R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
4. Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S. Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.
5. Taber, L.H.V. Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B. Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. J.Clin. Microbiol. 257 : 763-767.
7. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
8. Infektionen durch Adenovirus - Infektionskrankheiten - MSD Manual Profi-Ausgabe (msdmanuals.com), as of August 8th, 2024
9. Adenoviren: www.infektionsschutz.de, as of August 8th, 2024

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15. Lähteet

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
3. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
4. Hall, C.B.,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. J.Clin. Microbiol. 257 : 763-767.
7. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
8. Infektionen durch Adenovirus - Infektionskrankheiten - MSD Manual Profi-Ausgabe (msdmanuals.com), as of August 8th, 2024
9. Adenoviren: www.infektionsschutz.de, as of August 8th, 2024

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15. Referințe

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C. Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
3. Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
4. Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S. Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N. Engl. J. Med. 308 : 1443-1447.
5. Taber, L.H.V., Knight, B.E., Gilbert, H.W., McClung, S.Z., Wilson, H.J., Norton, J.M., Thurston, W.H., Gordon, R.L., Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
6. Ahluwalia, G.J., Embree, P., McNicol, B., Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. J. Clin. Microbiol. 257 : 763-767.
7. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
8. Infektionen durch Adenovirus - Infektionskrankheiten - MSD Manual Profi-Ausgabe (msdmanuals.com), as of August 8th, 2024
9. Adenoviren: www.infektionsschutz.de, as of August 8th, 2024

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**VivaChek Biotech (Hangzhou) Co., Ltd.**
Lotus NL B.V.