NADAL® Gonorrhea Test (test cassette)

REF 611905N-10

de	Gebrauchsanweisung	2	cs	Návod k použití	33
en	Instructions for use	7	nl	Gebruiksaanwijzing	37
fr	Instructions d'utilisation	11	bg	Инструкции за употреба	41
es	Instrucciones de uso	16		Symbols	47
it	Istruzioni per l'uso	20		Our Teams	48
pl	Sposób użycia	24			
pt	Instruções de Utilização	29			

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15. Referenzen

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC, 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6)
- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* und *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.

5. National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
6. Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Productive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
7. Ebert JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.
8. Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," Acta Cytol, 1985, 29:1043-6.
9. Gupta PK, "Microbiology, Inflammation, and Viral Infection," Comprehensive Cytopathology, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-52.
10. Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," Acta Cytol, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," Diagn Cytopathol, 1988, 4:169-76.
11. Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," J Clin Pathol, 1990, 43(9):758-61.
12. Schachter J, "Chlamydiae," Manuel of Clinical Laboratory Immunology, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

Rev.0, 2022-05-23 OM

1. Intended Use

The NADAL® Gonorrhea Test is a lateral flow chromatographic immunoassay for the qualitative detection of *Neisseria gonorrhoeae* antigens in female endocervical or male urethral swab specimens. The test is intended for use as an aid in the diagnosis of gonococcal infections (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® Gonorrhea Test is designed for professional use only.

2. Introduction and Clinical Significance

Gonorrhea is a sexually transmitted disease caused by the bacterium *Neisseria gonorrhoeae*. Gonorrhea is one of the most common infectious bacterial diseases and is most frequently transmitted during sexual intercourse, including vaginal, oral and anal sex. The pathogen can infect the throat, causing a severe sore throat. It can also infect the anus and rectum, triggering a condition called proctitis. In females, it can infect the vagina, causing irritation and discharge (vaginitis). Infection of the urethra may cause urethritis with burning, painful urination and discharge. If women have symptoms, they often notice vaginal discharge, increased urinary frequency and urinary discomfort. The spread of bacteria to the fallopian tubes and abdomen may cause severe lower-abdominal pain and fever. The average incubation period of gonorrhea is approximately 2 to 5 days following sexual contact with an infected partner. Nevertheless, symptoms may appear as late as 2 weeks afterwards. However, there are 5%-20% of men and 60% of women that do not show any symptoms. A preliminary diagnosis of gonorrhea can be made at the time of examination. In women, gonorrhea is a common cause of pelvic inflammatory disease (PID). PID can lead to internal abscesses and long-lasting, chronic pelvic pain. PID can damage the fallopian tubes enough to cause infertility or increase the risk of ectopic pregnancy.

3. Test Principle

The NADAL® Gonorrhea Test enables the detection of *Neisseria gonorrhoeae* antigens through the visual interpretation of colour development on the internal test strip. Specific anti-*Neisseria gonorrhoeae* antibodies are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with anti-*Neisseria gonorrhoeae* antibodies which are conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there is a sufficient number of *Neisseria gonorrhoeae* antigens in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® Gonorrhea test cassettes
- 10 extraction tubes, incl. dropper caps

Additional material provided in accordance with 93/42/EEC:

- 10 sterile swabs (female), CE 0197

Jiangsu Hanheng Medical Technology Co. Ltd.,
16-B4, #1, North Qingyang Road, Tianning District,
213017, Changzhou, Jiangsu, China (authorised
EU representative: Luxus Lebenswelt GmbH,
Kochstr. 1, 47877 Willich, Germany)

- 1 bottle of Reagent 1 (white cap), (4 mL)
- 1 bottle of Reagent 2 (red cap)*, (4 mL)
- 1 reagent holder
- 1 package insert

*contains the following preservative: ProClin™ 300: <0.09%.

5. Additional Materials Required

- Timer
- Sterile urethral swabs (male)
- Optional: Separate swab or cotton ball

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use Reagents 1 & 2 if they are discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.

- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Do not use swabs from damaged pouches.
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

Note: The quality of specimens obtained is of extreme importance. Detection of *Neisseria gonorrhoeae* requires a rigorous and thorough collection technique which provides cellular material rather than just body fluids. Do not use 0.9% sodium chloride to treat swabs before collecting specimens.

Female endocervical specimen:

Prior to specimen collection, remove excess mucus from the endocervical area using a separate swab or a cotton ball (not provided with the test kit). Dispose of the swab or the cotton ball.

Insert the provided swab into the endocervical canal past the squamocolumnar junction. This enables to collect columnar or cuboidal epithelial cells, which are the main reservoir of *Neisseria gonorrhoeae*.

Firmly rotate the swab for 15-20 seconds, avoiding contamination with exocervical or vaginal cells.

If the swab is to be tested immediately, place it into a tube.

Male urethral specimen:

Note: Instruct the patients not to urinate for at least 2 hours prior to specimen collection.

Use a standard wire-shafted fiber-tipped swab for urethral specimen collection.

Insert the swab 2-4 cm into the urethra, rotate for 3-5 seconds and remove it.

If the swab is to be tested immediately, place the swab into a tube.

Note: Do not place a swab into a tube containing transport medium. Transport medium may interfere with the assay. Viability of bacteria is not required for the assay. Swab specimens should be tested immediately after collection. For best test performance, use freshly collected specimens. If a swab cannot be processed immediately, it should be placed into a dry transport tube for storage or transportation. Swabs can be stored at room temperature (15-30°C) for up to 4-6 hours or at 2-8°C for up to 24 hours. Do not freeze specimens. All specimens should be brought to room temperature (15-30°C) prior to testing.

Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of specimens and lead to inaccurate test results.

To ensure optimal performance, use the female swabs supplied in the test kit only.

9. Test Procedure

Bring tests, specimens, reagents and/or controls to room temperature (15-30°C) prior to testing.

To avoid cross contamination, do not allow the tips of the reagent bottles to come into contact with specimen material.

1. Place a clean extraction tube labeled with the patient or control identification into the designated area of the reagent holder.
2. Holding the bottle vertically, add 6 drops of reagent 1 to the extraction tube. Then close the bottle.
3. Insert the swab into the extraction tube. Swirl the swab and leave it in the tube for 2 minutes.
4. Holding the bottle vertically, add 8 drops of reagent 2 to the extraction tube and mix well. Then close the bottle.
5. Swirl the swab for about 1 minute.
6. Remove the swab, pressing it firmly against the wall of the tube and compressing the wall of the tube against the swab to expunge as much liquid as possible from the swab.
7. Dispose of the swab in accordance with guidelines for the handling of infectious agents.
8. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
9. Place the test cassette on a clean and level surface.

10. Attach a dropper cap to the tube.
11. Invert the tube and transfer 3 drops (approximately 110 µL) of the extracted solution from the extraction tube to the specimen well (S) of the test cassette.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

12. Start the timer.


As the test begins to run, you will observe a coloured liquid migrate along the membrane.

13. Wait for the coloured line(s) to appear. Read the test result after 15 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No coloured line develops in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test cassette: A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Gonorrhea Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of *Neisseria gonorrhoeae* antigens in female endocervical or male urethral swab specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of *Neisseria gonorrhoeae* antigens can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- Detection of *Neisseria gonorrhoeae* depends on the number of bacteria present in the specimen. It may be affected by patient factors, such as age, a history of STDs, the presence of symptoms, etc. The lower detection limit of this test may vary according to the serovar.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods (e.g. bacterial culture or nucleic acid amplification tests (NAATs)) is recommended. A negative result does not at any time preclude the possibility of a gonococcal infection.
- The NADAL® Gonorrhea Test only detects the presence of antigens from both viable and non-viable *Neisseria gonorrhoeae* in specimens and should not be used as the sole criterion for a diagnosis of a gonococcal infection.

- Positive results do not preclude co-infections with other pathogens.
- Female patients with vaginal discharge should be assessed for risk factors of cervicitis and pelvic inflammatory disease caused by other microorganisms including *Candida albicans*, *Trichomonas vaginalis* or bacterial vaginosis.
- Therapeutic success or failure cannot be determined using the test, as *Neisseria gonorrhoeae* antigens may persist following appropriate antimicrobial therapy.
- Excessive blood (>50 µL in case of female endocervical swab specimens and >20 µL in case of male urethral swab specimens) may cause false positive results. Endocervical swab specimens from female patients should not be collected during menstrual period.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® Gonorrhea Test was evaluated using clinical swab specimens whose status was confirmed using bacterial culture.

The results are presented in the following tables:

Female endocervical swab specimens:

NADAL® Gonorrhea Test	Bacterial culture			
		Positive	Negative	Total
	Positive	62	21	83
	Negative	2	473	475
Total		64	494	558

Diagnostic sensitivity: 96.9% (89.2% - 99.6%)*

Diagnostic specificity: 95.8% (93.6% - 97.4%)*

Overall agreement: 95.9% (93.9% - 97.4%)*

*95% confidence interval

Male urethral swab specimens:

NADAL® Gonorrhea Test	Bacterial culture			
		Positive	Negative	Total
	Positive	134	2	136
	Negative	3	389	392
Total		137	391	528

Diagnostic sensitivity: 97.8% (93.7% - 99.6%)*

Diagnostic specificity: 99.5% (98.2% - 99.9%)*

Overall agreement: 99.1% (97.8% - 99.7%)*

*95% confidence interval

Analytical performance

Analytical specificity

Cross-reactivity study

Buffer spiked with the following potentially cross-reacting microorganisms at the concentration of 1.0×10^7 organisms/mL showed no cross-reactivity when tested using the NADAL® Gonorrhea Test:

<i>Acinetobacter calcoaceticus</i>	<i>Pseudomonas aeruginosa</i>
<i>Acinetobacter spp.</i>	<i>Gardnerella vaginalis</i>
<i>Enterococcus faecalis</i>	<i>Salmonella enterica</i>
<i>Enterococcus faecium</i>	<i>Proteus vulgaris</i>
<i>Staphylococcus aureus</i>	<i>Haemophilus influenzae</i>
<i>Proteus mirabilis</i>	<i>Klebsiella pneumoniae</i>
<i>Chlamydia trachomatis</i>	<i>Ureaplasma urealyticum</i>
<i>Streptococcus, group B/C</i>	<i>Trichomonas vaginalis</i>
<i>Mycoplasma hominis</i>	<i>Candida albicans</i>

Interference study

Negative controls, as well as positive controls prepared by diluting *Neisseria gonorrhoeae* culture with negative buffer were spiked with the following potentially interfering substances, showed no interference with the NADAL® Gonorrhea Test:

Paracetamol (Acetaminophen)	10 mg/mL
Aspirin (acetylsalicylic acid)	20 mg/mL
Ibuprofen	10 mg/mL
Acidic urine (pH 4)	/
Alkaline urine (pH 9)	/
Leukocytes	1 x 10 ⁶ cells/mL
Cellular debris	/
Vitamin C	20 µg/mL
Calciumphosphat	1.5% W/V
Albumin	0.2 mg/mL
Bilirubin	0.2 mg/mL
Glucose	2 µg/mL
Nitrate	500 mg/L
Ketones	1.6 mg/mL
Urobilinogen	20 mg/mL
10% blood	/
30% blood	/
Contraceptive gel	0.25% W/V
Spermicide	0.25% W/V
Moisturiser	0.25% W/V
Haemorrhoidal anaesthetic	0.25% W/V
Body oil	1 mg/mL
Powder	10 mg/mL
Antifungal cream	0.25% W/V
Vaginal lubricants	0.25% W/V
Feminine spray	0.25% W/V

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative, low, medium and high positive controls using 3 lots of the NADAL® Gonorrhea tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5% - 100% for each lot). The NADAL® Gonorrhea Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of negative, low, medium and high positive controls. Testing was performed by 3 operators using 3 independent NADAL® Gonorrhea test lots at 3 different sites on 5 separate days. >99% of the specimens were correctly identified (225/225 correct tests per concentration, 95% confidence interval: 99.6% - 100%). The NADAL® Gonorrhea Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® Gonorrhea Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC, 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6).
- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* and *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.
- National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
- Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Reproductive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
- Ehert JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.
- Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," Acta Cytol, 1985, 29:1043-6.
- Gupta PK, "Microbiology, Inflammation, and Viral Infection," Comprehensive Cytopathology, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-52.
- Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," Acta Cytol, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," Diagn Cytopathol, 1988, 4:169-76.
- Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," J Clin Pathol, 1990, 43(9):758-61.
- Schachter J, "Chlamydiae," Manuel of Clinical Laboratory Immunology, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

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15. Bibliographie

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC, 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6)
- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* et *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.
- National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
- Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Reproductive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
- Ehrert JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.
- Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," Acta Cytol, 1985, 29:1043-6.

9. Gupta PK, "Microbiology, Inflammation, and Viral Infection," Comprehensive Cytopathology, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-52.
10. Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," Acta Cytol, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," Diagn Cytopathol, 1988, 4:169-76.
11. Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," J Clin Pathol, 1990, 43(9):758-61.
12. Schachter J, "Chlamydiae," Manuel of Clinical Laboratory Immunology, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

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Jiangsu Hanheng Medical Technology Co. Ltd.,
16-B4, #1, North Qingyang Road, Tiansheng District,
213017, Changzhou, Jiangsu, China (representante EU autorizado: Luxus Lebenswelt GmbH, Kochstr. 1, 47877 Willich, Germany)

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15. Referencias

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC, 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6)
- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* and *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.
- National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
- Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Reproductive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
- Ehrt JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.
- Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," Acta Cytol, 1985, 29:1043-6.
- Gupta PK, "Microbiology, Inflammation, and Viral Infection," Comprehensive Cytopathology, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-52.
- Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," Acta Cytol, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," Diagn Cytopathol, 1988, 4:169-76.
- Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," J Clin Pathol, 1990, 43(9):758-61.
- Schachter J, "Chlamydiae," Manual of Clinical Laboratory Immunology, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

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Jiangsu Hanheng Medical Technology Co. Ltd.,
16-B4, #1, North Qingyang Road, Tianning District,
213017, Changzhou, Jiangsu, Cina (rappre-
sentante autorizzato UE: Luxus Lebenswelt GmbH,
Kochstr. 1, 47877 Willich, Germany)

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15. Bibliografia

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC, 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(1RR-6).
- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* and *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.
- National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
- Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Reproductive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
- Ehert JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.
- Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," Acta Cytol, 1985, 29:1043-6.
- Gupta PK, "Microbiology, Inflammation, and Viral Infection," Comprehensive Cytopathology, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-52.
- Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," Acta Cytol, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," Diagn Cytopathol, 1988, 4:169-76.
- Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," J Clin Pathol, 1990, 43(9):758-61.
- Schachter J, "Chlamydiae," Manuel of Clinical Laboratory Immunology, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

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15. Bibliografia

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC, 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6).
- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* and *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.
- National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
- Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Reproductive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
- Ehert JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.

8. Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," *Acta Cytol*, 1985, 29:1043-6.
9. Gupta PK, "Microbiology, Inflammation, and Viral Infection," *Comprehensive Cytopathology*, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-52.
10. Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," *Acta Cytol*, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," *Diagn Cytopathol*, 1988, 4:169-76.
11. Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," *J Clin Pathol*, 1990, 43(9):758-61.
12. Schachter J, "Chlamydiae," *Manuel of Clinical Laboratory Immunology*, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

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15. Reference

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC., 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6)

- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* and *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.
- National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
- Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Reproductive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
- Ehert JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.
- Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," Acta Cytol, 1985, 29:1043-6.
- Gupta PK, "Microbiology, Inflammation, and Viral Infection," Comprehensive Cytopathology, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-52.
- Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," Acta Cytol, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," Diagn Cytopathol, 1988, 4:169-76.
- Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," J Clin Pathol, 1990, 43(9):758-61.
- Schachter J, "Chlamydiae," Manuel of Clinical Laboratory Immunology, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

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15. Referenties

- Knapp, J.S. et al. *Neisseria gonorrhoeae*. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC, 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6).
- Forbes B.A., Sahm D.F., Weissfeld A.S. *Neisseria* and *Moraxella catarrhalis*. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.
- National Institute of Allergy and Infectious Diseases, National Institute of Health, US Department of Health and Human Services, NIAID Fact Sheet on Gonorrhea, October 2004.
- Bibbo M and Wied GL, "Inflammation Reaction and Microbiology of the Female Reproductive Tract," Compendium on Diagnostic Cytology, Wied GL, Keebler CM, Koss LG, et al, eds, Chicago, IL: Tutorials of Cytology, 1992, 63-8.
- Ehert JM and Judson FN, "Genital Chlamydia Infections," Clin Lab Med, 1989, 9(3):481-500.
- Gay JD, Donaldson LD, and Goellner JR, "False-Negative Results in Cervical Cytologic Studies," Acta Cytol, 1985, 29:1043-6.
- Gupta PK, "Microbiology, Inflammation, and Viral Infection," Comprehensive Cytopathology, 1st ed, Bibbo M, ed, Philadelphia, PA: WB Saunders Co, 1991, 115-152.
- Luff RD, "The Bethesda System for Reporting Cervical/Vaginal Diagnosis," Acta Cytol, 1993, 37(2): 115-24. Maguire NC, "Current Use of the Papanicolaou Class System in Gynecologic Cytology," Diagn Cytopathol, 1988, 4:169-76.
- Miettinen A, Heinonen PK, Teisala K, et al, "Antigen-Specific Serum Antibody Response to Chlamydia trachomatis in Patients With Pelvic Inflammatory Disease," J Clin Pathol, 1990, 43(9):758-61.
- Schachter J, "Chlamydiae," Manuel of Clinical Laboratory Immunology, 4th ed, Vol 2, Chapter 96, Rose NR, Conway de Macario E, Fahey JL, et al, eds, Washington, DC: American Society for Microbiology, 1992, 661-6.

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- Knapp, J.S. et al. Neisseria gonorrhoeae. Manual of Clinical Microbiology, Sixth Edition, ASM Press, Washington DC., 324-325 (1995).
- Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines 2002. Morbidity and Mortality Weekly Report (2002), 51(RR-6)
- Forbes B.A., Sahm D.F., Weissfeld A.S. Neisseria and Moraxella catarrhalis. Bailey & Scott's Diagnostic Microbiology, Tenth Edition, Mosby, St. Louis, 597-605 (1998).
- Summary of the Notifiable Diseases, United States, 1998, Morbidity and Mortality Weekly Report (1999), 47(53): 1-93.

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