NADAL® HAV IgG/IgM Test (test cassette)

REF 622970N-30

de	Gebrauchsanweisung	2	pt	Instruções de Utilização	26
en	Instructions for use	6	cs	Návod k použití	30
fr	Instructions d'utilisation	10		Symbols	35
es	Instrucciones de uso	14		Our Teams	36
it	Istruzioni per l'uso	18			
pl	Sposób użycia	22			

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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15. Referenzen

1. Feinstone SM, Kapikian AZ, Purcell RH. Hepatitis A: Detection by immune electron microscopy of a viruslike antigen associated with acute illness. Science. 1973;182:1026-1028.
2. Matheny SC, Kingery JE. Hepatitis A. Am Fam Physician. 2012 Dec 1;86(11):1027-34.
3. Rose NR, Hamilton RG, Detrick B. Manual of Clinical Laboratory Immunology. American Society for Microbiology, Washington DC, 2002
4. Nainan OV, Xia G, Vaughan G, Margolis HS. Diagnosis of hepatitis A virus infection: a molecular approach. Clin Microbiol Rev. 2006 Jan;19(1):63-79.
5. Stapleton JT. Host immune response to hepatitis A virus. J Infect Dis, 1995, 171(Suppl 1):S9-S14.

Rev. 0, 2022-05-23 OM

1. Intended Use

The NADAL® HAV IgG/IgM Test is a lateral flow chromatographic immunoassay for the qualitative detection and differentiation of anti-hepatitis A virus IgG and IgM (anti-HAV IgG and IgM) in human whole blood, serum or plasma specimens. The test is intended for use as an aid in the diagnosis of active and/or past HAV infections and provides an initial presumptive diagnosis (see section 12 'Limitations'). Positive results must be confirmed using (an) alternative testing method(s) (such as PCR or ELISA) and clinical findings. The test procedure is not automated and requires no special training or qualification. The NADAL® HAV IgG/IgM Test is designed for professional use only.

2. Introduction and Clinical Significance

First isolated in 1973,¹ the hepatitis A virus (HAV) is a non-enveloped RNA virus of the family *Picornaviridae*, genus *Hepatovirus*. Its transmission occurs primarily through serial transmissions from person to person by the fecal-oral route. However, HAV can also be contracted from contaminated water, close personal contact, sexual contact, and illicit drug use. Following an incubation period of approximately 28 days (15-50 days), the disease is characterised by a rapid onset of symptoms. Signs and symptoms can include nausea, vomiting, diarrhoea, dark urine, jaundice, fever, headaches, weight loss and abdominal pain. The likelihood of symptoms increases with age. Normally, acute illness does not last more than two months. There is no chronic viral shedding and no chronic stage of the disease, though recurrences, acute fulminant hepatitis and other complications may occur.² Anti-HAV IgM are detectable at or prior to the onset of clinical illness and decline in about 3 to 6 months. Anti-HAV IgG appear soon after IgM, persist for years after infection and confer lifelong immunity.³ The presence of anti-HAV IgG and the absence of anti-HAV IgM can be used to differentiate between past and current infections. IgM assays can detect antibodies in individuals who have recently received a hepatitis A vaccine for a short period of time. However, lower concentrations detected 4 to 6 months after the onset of infection typically do not produce a positive test result.^{4,5}

3. Test Principle

The NADAL® HAV IgG/IgM Test is a lateral flow chromatographic immunoassay for the qualitative detection and differentiation of anti-HAV IgG and IgM in human whole blood, serum or plasma specimens.

Mouse anti-human IgG and mouse anti-human IgM are immobilised in the test line regions of the membrane labelled 'IgG' and 'IgM' respectively.

During the test, anti-HAV IgG and/or anti-HAV IgM, if present in the specimen, bind to the HAV antigens conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. The complexes are then captured by the mouse anti-human IgG and mouse anti-human IgM in the test line regions 'IgG' and 'IgM' respectively. The presence of a coloured line in the test line region 'IgG' and/or 'IgM' indicates a positive anti-HAV IgG and/or anti-HAV IgM result. The absence of a coloured line in the test line region 'IgG' and/or

'IgM' indicates a negative anti-HAV IgG and/or anti-HAV IgM result.

The formation of a coloured line in the control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 30 NADAL® HAV IgG/IgM test cassettes
- 30 capillary tubes (5 µL)
- 1 buffer (6 mL)*
- 1 package insert

*containing the following preservative: ProClin™ 300: 0.05%.

No hazard labelling for ProClin™ 300 is required according to Regulation (EC) N° 1272/2008 CLP.

5. Additional Materials Required

- Specimen collection containers (appropriate for specimen material to be tested)
- Centrifuge (for serum or plasma specimens only)
- Alcohol pads
- Lancets (for fingerstick whole blood specimens only)
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Do not use the test if the foil pouch does not contain a desiccant bag or if a desiccant bag is damaged inside the foil pouch.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.

- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

The NADAL® HAV IgG/IgM Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

Venipuncture whole blood specimens

Containers containing anticoagulants, such as EDTA, sodium citrate or heparin sodium should be used for the preparation of venous whole blood or plasma specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time.

If the test is to be run within 24 hours of specimen collection, whole blood collected by venipuncture should be stored at 2-8°C.

Do not freeze whole blood specimens.

Serum and plasma specimens

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 5 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well as

well as carefully prior to testing. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Icteric, lipemic, haemolysed, heat-treated and contaminated specimens may lead to inaccurate test results.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

2. Place the test cassette on a clean and level surface.

3. Holding the capillary tube vertically, draw the specimen (whole blood/ serum/ plasma) up to the fill line (approximately 5 µL) and add it to the to the specimen well (S) of the test cassette.

Alternatively, a micropipette (5 µL) may be used.


Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

4. Holding the buffer bottle vertically, add 2 drops (approximately 60-80 µL) of buffer to the buffer well (B).

Avoid trapping air bubbles in the buffer well (B) and do not add any solution to the result area.

5. Start the timer.


6. Wait for the coloured line(s) to appear. Read the test result after 15 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation


Positive for IgM:

A coloured line develops in the control line region 'C'. Another coloured line develops in the test line region 'IgM'. The test result indicates the presence of anti-HAV IgM.

Positive for IgG:

A coloured line develops in the control line region 'C'. Another coloured line develops in the test line region 'IgG'. The test result indicates the presence of anti-HAV IgG.

Positive for IgM and IgG:

A coloured line develops in the control line region 'C'. A coloured line develops in the test line region 'IgM' and another one develops in the test line region 'IgG'. The

test result indicates the presence of anti-HAV IgM and IgG.

Note: The colour intensity in the test line regions 'IgG' or 'IgM' may vary depending on the concentration of anti-HAV antibodies in the specimen. Therefore, any shade of colour in the test line region 'IgG' or 'IgM' should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region 'C'. No coloured lines develop in the test line regions 'IgM' or 'IgG'.

Invalid:

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region 'C' is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® HAV IgG/IgM Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection and differentiation of anti-HAV IgG and IgM in human whole blood, serum or plasma specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of anti-HAV antibodies can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factor (≥ 1000 IU/mL) may affect test results.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.
- The NADAL® HAV IgG/IgM Test only detects the presence of anti-HAV antibodies in specimens and should not be used as the sole criterion for a diagnosis of an HAV infection.

- An HAV infection may progress rapidly. If the test result is negative but clinical symptoms persist, additional testing using other clinical methods such as PCR or ELISA is recommended.
- A negative result does not at any time preclude the possibility of exposure to or an acute infection with HAV. Specimens collected too early in the course of an infection may contain IgM levels which are below the detection limit of this test.
- A negative test result may occur if the titer of HAV antibodies present in the specimen is below the detection limit of the test or if the detected HAV antibodies were not present during the stage of disease in which the specimen was collected.
- The presence of detectable anti-HAV IgG may be due to past infections and/or vaccination.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® HAV IgG/IgM Test was evaluated using clinical plasma specimens whose status was confirmed using ELISA.

The results are presented in the following tables:

NADAL® HAV IgG/IgM Test (IgG) vs. ELISA

NADAL® HAV IgG/IgM Test (IgG)	anti-HAV IgG ELISA		
	Positive	Negative	Total
Positive	125	4	129
Negative	0	71	71
Total	125	75	200

Diagnostic sensitivity: >99.9% (97.1% - 100.0%)*

Diagnostic specificity: 94.7% (86.9% - 98.5%)*

Overall agreement: 98.0% (95.0% - 99.5%)*

*95% confidence interval

NADAL® HAV IgG/IgM Test (IgM) vs. ELISA

NADAL® HAV IgG/IgM Test (IgM)	anti-HAV IgM ELISA		
	Positive	Negative	Total
Positive	91	7	98
Negative	5	203	208
Total	96	210	306

Diagnostic sensitivity: 94.8% (88.3% - 98.3%)*

Diagnostic specificity: 96.7% (93.3% - 98.7%)*

Overall agreement: 96.1% (93.3% - 98.0%)*

*95% confidence interval

Performance on Boston Biomedica Inc (BBI) Anti-HAV Seroconversion Panel

The performance of the NADAL® HAV IgG/IgM Test was evaluated using the BBI Anti-HAV Seroconversion Panel (PHT903). The results are presented in the following table:

BBI Reference Panel: HAV ELISA		NADAL® HAV IgG/IgM Test	
Type	Number	Agreement	
HAV IgG positive	8	8 (100%)	
HAV IgG negative	2	2 (100%)	
HAV IgM positive	8	8 (100%)	
HAV IgM negative	2	2 (100%)	

Analytical performance

Detection limit

The detection limit of the NADAL® HAV IgG/IgM Test is 70 mIU/mL for IgG (NIBSC 97/646).

Analytical specificity

Cross-reactivity

Serum specimens of the following infectious diseases and immunomarkers (anti-HAV IgG as well as anti-HAV IgM negative) were tested using the NADAL® HAV IgG/IgM Test.

HBV	Malaria
HCV	Tuberculosis
HEV	<i>T. pallidum</i>
HIV	Typhoid fever
hCG	ANA
Dengue	HAMA
<i>H. pylori</i>	Rheumatoid factors (up to 1000 IU/mL)

No cross-reactivity with the specimens was observed when tested using the NADAL® HAV IgG/IgM Test.

Interfering substances

Negative and low positive anti-HAV IgG as well as anti-HAV IgM specimens spiked with the following potentially interfering substances showed no interference with the NADAL® HAV IgG/IgM Test.

Substance	Concentration
Albumin	60 g/L
Bilirubin	0.2 mg/mL
Creatinine	442 µmol/L
EDTA	3.4 µmol/L
Glucose	55 mmol/L
Haemoglobin	2 g/L
Heparin sodium	3000 IU/mL
Salicylic acid	4.24 mmol/L
Sodium citrate	3.4%

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative, low, medium and high positive anti-HAV IgG and low positive anti-HAV IgM serum specimens. Testing was performed by one operator using 3 lots of the NADAL® HAV IgG/IgM tests at one site on one day. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 96.4% - 100% for each lot). The NADAL® HAV IgG/IgM Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 2 replicates of negative, low, medium and high positive anti-HAV IgG and low positive anti-HAV IgM serum specimens. Testing was performed by 2 operators using 3 independent NADAL® HAV IgG/IgM test lots at 3 different sites on 5 separate days. >99% of the specimens were correctly identified (180/180 correct tests per concentration, 95% confidence interval: 99.6% - 100%). The NADAL® HAV IgG/IgM Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® HAV IgG/IgM Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Feinstone SM, Kapikian AZ, Purcell RH. Hepatitis A: Detection by immune electron microscopy of a viruslike antigen associated with acute illness. *Science*. 1973;182:1026-1028.
2. Matheny SC, Kingery JE. Hepatitis A. *Am Fam Physician*. 2012 Dec 1;86(11):1027-34.
3. Rose NR, Hamilton RG, Detrick B. *Manual of Clinical Laboratory Immunology*. American Society for Microbiology, Washington DC, 2002
4. Nainan OV, Xia G, Vaughan G, Margolis HS. Diagnosis of hepatitis a virus infection: a molecular approach. *Clin Microbiol Rev*. 2006 Jan;19(1):63-79.
5. Stapleton JT. Host immune response to hepatitis A virus. *J Infect Dis*. 1995; 171(Suppl 1):S9-S14.

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15. Bibliographie

- Feinstone SM, Kapikian AZ, Purcell RH. Hepatitis A: Detection by immune electron microscopy of a viruslike antigen associated with acute illness. *Science*. 1973;182:1026-1028.
- Matheny SC, Kingery JE. Hepatitis A. *Am Fam Physician*. 2012 Dec 1;86(11):1027-34.
- Rose NR, Hamilton RG, Detrick B. *Manual of Clinical Laboratory Immunology*. American Society for Microbiology, Washington DC, 2002
- Nainan OV, Xia G, Vaughan G, Margolis HS. Diagnosis of hepatitis A virus infection: a molecular approach. *Clin Microbiol Rev*. 2006 Jan;19(1):63-79.
- Stapleton JT. Host immune response to hepatitis A virus. *J Infect Dis*, 1995, 171(Suppl 1):S9-S14.

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15. Referencias

- Feinstone SM, Kapikian AZ, Purcell RH. Hepatitis A: Detection by immune electron microscopy of a viruslike antigen associated with acute illness. *Science*. 1973;182:1026-1028.
- Matheny SC, Kingery JE. Hepatitis A. *Am Fam Physician*. 2012 Dec 1;86(11):1027-34.
- Rose NR, Hamilton RG, Detrick B. *Manual of Clinical Laboratory Immunology*. American Society for Microbiology, Washington DC, 2002
- Nainan OV, Xia G, Vaughan G, Margolis HS. Diagnosis of hepatitis a virus infection: a molecular approach. *Clin Microbiol Rev*. 2006 Jan;19(1):63-79.
- Stapleton JT. Host immune response to hepatitis A virus. *J Infect Dis*. 1995, 171(Suppl 1):S9-S14.

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15. Bibliografia

1. Feinstone SM, Kapikian AZ, Purcell RH. Hepatitis A: Detection by immune electron microscopy of a viruslike antigen associated with acute illness. Science. 1973;182:1026-1028.
2. Matheny SC, Kingery JE. Hepatitis A. Am Fam Physician. 2012 Dec 1;86(11):1027-34.
3. Rose NR, Hamilton RG, Detrick B. Manual of Clinical Laboratory Immunology. American Society for Microbiology, Washington DC, 2002
4. Nainan OV, Xia G, Vaughan G, Margolis HS. Diagnosis of hepatitis a virus infection: a molecular approach. Clin Microbiol Rev. 2006 Jan;19(1):63-79.
5. Stapleton JT. Host immune response to hepatitis A virus. J Infect Dis, 1995, 171(Suppl 1):S9-S14.

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15. Referências

- Feinstone SM, Kapikian AZ, Purcell RH. Hepatitis A: Detection by immune electron microscopy of a viruslike antigen associated with acute illness. Science. 1973;182:1026-1028.
- Matheny SC, Kingery JE. Hepatitis A. Am Fam Physician. 2012 Dec 1;86(11):1027-34.
- Rose NR, Hamilton RG, Detrick B. Manual of Clinical Laboratory Immunology. American Society for Microbiology, Washington DC, 2002
- Nainan OV, Xia G, Vaughan G, Margolis HS. Diagnosis of hepatitis A virus infection: a molecular approach. Clin Microbiol Rev. 2006 Jan;19(1):63-79.
- Stapleton JT. Host Immune response to hepatitis A virus. J Infect Dis, 1995, 171(Suppl 1):S9-S14.

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